

## Now you can keep on going

You have already seen and done a lot. But your journey of discovery isn't over yet. Effortless staying power. Live life to the full. That's what Rollz is all about. Which is why we have designed the modern Rollz Motion, a rollator that converts into a wheelchair at a flick of the wrist. And of course, just as easily back again. Rollz takes you further.

Be it for an outdoor activity, a long holiday, a treasured exposition or a shopping afternoon, this 2-in-1 rollator inspires you to keep on going without effort. In addition the Rollz Motion folds in a compact package and fits into any car.

Rollz  
takes you  
further



### Rollz International

Pedro de Medinalaan 63

1086 XP Amsterdam

The Netherlands

info@rollz.com

Facebook: Rollz International

Twitter: www.twitter.com/rollzint

[www.rollz.com](http://www.rollz.com)



Rollator and  
Wheelchair  
in one



2012



reddot design award  
honourable mention 2012



# Rollator and Wheelchair

The Rollz Motion is an ergonomically designed rollator with a strong but lightweight frame which enables you to walk steadily, upright and comfortably. And if you want to rest along the way, the Rollz Motion converts into a wheelchair at a flick of the wrist. And of course, just as easily back again.



## From rollator to wheelchair at a flick of the wrist

Before you leave home, click the Rollz wheelchair attachment into place on the rollator.

1. Pull the backrest into position in one flowing motion.
2. Turn the handles from the rollator position into the push position, and adjust them to the right height for the person who's going to push.
3. Then click the footrests into place at the front.
4. And continue your journey in comfort and style.

## Product Details

### Compact

The Rollz Motion is easily foldable into a compact package.



### Ease

Cross doorsteps and curbs easily with the threshold aid.



### Stable

Walking tall, steadily and supported between the ergonomic hand grips that are adjustable in height.



### Safe

Innovative drum brakes for safe and balanced braking.

### Convenient

Resting on a comfortable cushioned seat.

### Comfortable

Comfortable ride and natural suspension thanks to the foam-filled rear tyres.

### Manoeuvrable

Light steering and easy motion with the manoeuvrable front wheels.

## Facts & Figures Rollz Motion

### Weight

Rollz Motion: 10.6 kg

Rollz wheelchair package: 3.7 kg

### User

Maximum weight: 125 kg

Minimum and maximum height: 1.60m to 1.90m

### Dimensions

Rollz Motion:

Unfolded HxWxL cm: 86-102 x 65.5 x 63

Folded HxWxL cm: 86-102 x 29 x 63

Wheelchair seat:

HxWxL cm: 55 x 45 x 38.5

### Accessories

Rollz Shopper

Rollz cup holder

Rollz cane holder

Rollz Motion is available in four colours:

orange, ice-blue, sand, cool grey



For information on Rollz dealers visit [www.rollz.com](http://www.rollz.com)

